**EmCyte Corporation**  
4331 Veronica S. Shoemaker Blvd., Suite 4  
Fort Myers, FL 33916

PL-42-13  
REV 11  
Date: August 2017

**PureBMC®**  
**Platelet Concentrating System 60mL**  
**PN: BC60-PURE**  
**Instruction for Use**

**ATTENTION OPERATING SURGEON**

**NOTE: DEVICE IS FOR SINGLE USE ONLY.** Discard the entire disposable system after one use, using an acceptable disposal method for products potentially contaminated with blood.

**DESCRIPTION**

PureBMC® is manufactured by EmCyte Corporation and is part of the GenesisCS Component Concentrating System's & family of kits. The kit concentrates bone marrow aspirate components to form platelet rich plasma by density through the use of the concentrating device, concentrating accessory and the Executive Series Centrifuge II provided by EmCyte Corporation. The system contains syringes, a bone marrow needle accessory and the concentrating device accessories. The PureBMC® concentrating system permits the rapid preparation of platelet rich plasma from the patient's bone marrow aspirate drawn at the time of treatment.

**MATERIALS**

The materials used are syringes, needles, tubing, connectors, and concentrating devices. The materials consist of medical grade polymers, elastomers and stainless steel that are suitable for use in medical devices. All components in this system are packaged, labeled and sterilized as indicated by the manufacturer's labeling. All components in this system are latex-free.

**INDICATIONS FOR USE STATEMENTS**

**FDA 510(k) BK050055**

The GenesisCS Component Concentrating System is designed to be used for the safe and rapid preparation of autologous platelet rich plasma (PRP) from a small sample of blood at the patient's point of care. The PRP can be mixed with autograft and allograft bone prior to application to an orthopedic surgical site as deemed necessary by the clinical use requirements.

**FDA 510(k) K07066**

The GenesisCS Component Concentrating System is designed to be used in a clinical laboratory or intraoperatively at the point of care for the safe and rapid preparation of platelet poor plasma and platelet concentrate from a small sample of blood and for a preparation of a cell concentrate from bone marrow. The safety and effectiveness of this device for in vivo indications for use has not been established.

**FDA 510(k) K024080**

The SECQUIRE cell separator is designed for the safe and effective preparation of low volume platelet rich plasma, and platelet poor plasma at the point of care.

**CE0653**

The Concentrating System is designed to be used for the safe and rapid preparation of autologous platelet rich plasma (PRP) from a small sample of blood and for a preparation of a cell concentrate from bone marrow aspirate. They can be mixed with autograft and allograft bone prior to application to an orthopedic surgical site and can be used for wound care and other soft tissue applications as deemed necessary by the clinical use requirements.

**WARNING AND PRECAUTIONS**

1. Use proper safety precautions to guard against needle sticks.
2. Follow manufacturer instructions when using centrifuge. Use only EmCyte provided general purpose centrifuge. Outcomes using centrifuges from other manufacturers are unknown.
3. When using the bone marrow aspiration needle, follow manufacturer's instruction for use.
4. Do not use sterile components of this system if package is opened or damaged.
5. Single use device. Do not reuse.
6. Do not use after expiration date.
7. The surgeon is to be thoroughly familiar with the equipment and the surgical procedure prior to using this device.
8. Use prepared BMC, PPP or PRP within 4 hours after drawing blood and marrow from the patient, according to current AABB guidelines.
9. The patient is to be made aware of the general risks associated with bone marrow aspiration. These risks include, but are not limited to: hemorrhage, seroma formation, infection, and/or persistent pain at the site of aspiration.
10. The safety and effectiveness of this device for in vivo indications for use, such as bone healing and hemostasis, have not been established.
11. The BMC or PRP prepared by this device has not been evaluated for any clinical indications.
12. BMC prepared from bone marrow may contain higher levels of plasma free hemoglobin than PRP prepared from whole blood.
13. The BMC or PRP prepared by this device is NOT indicated for delivery to the patient's circulatory system.
14. Device is single use only. Do not attempt to clean or re-sterilize this product. After use, this product may be a potential biohazard.

**POSSIBLE ADVERSE EFFECTS**

1. Damage to blood vessels, hematoma, delayed wound healing and/or infection is associated with blood draw, bone marrow harvest, and/or surgical procedure.
2. Temporary or permanent nerve damage that may result in pain or numbness is associated with blood draw, bone marrow harvest, and/or surgical procedure.
3. Early or late postoperative infection is associated with surgical procedure.
4. Pain associated with site of bone marrow harvest.

**STERILITY**

The GenesisCS Component Concentrating System kits are sterilized by ETO exposure. Do not use any component from an opened or damaged package. Do not resterilize. Single Use Only. Do not use past the expiration date.

**INSTRUCTIONS FOR USE**

**NOTE:** Use standard aseptic technique throughout the following procedure.

**PureBMC®**  
**Platelet Concentrating System 60mL**

1. **BONE MARROW ASPIRATE DRAW:**  
Attach the sterile filter needle onto the 60mL syringe. Draw 10mL of Citrate Anticoagulant into the 60mL syringe. Remove the filter needle from the syringe. Attach the syringe to the BMA needle. Slowly draw 50mL of bone marrow aspirate from the patient filling the syringe to 60mL. Follow the BMA needle manufacturer's package insert to obtain bone marrow aspirate. Gently, but thoroughly mix the BMA and anticoagulant upon collection to prevent coagulation.
2. **FILTER:** Using the filter accessory, inject the BMA through the filter into a sterile 60mL syringe. Filter in the direction indicated on the filter.
3. **LOAD:** Slowly add the filtered and anticoagulated BMA into the Concentrating Device. Load into the Executive Series Centrifuge II machine.
4. **BALANCE:** Counterbalance the Concentrating Device. Place the counterbalance and Concentrating Device directly opposite to each other in the centrifuge rotor and make sure they contain the same amount of volume.
5. **FIRST SPIN:** Close centrifuge lid. Set RPM to 3800 and the time to 2.5 minutes. Press the start button. Once the centrifuge stops, remove the Concentrating Device.
6. **FIRST EXTRACTION & TRANSFER:** Attach the Aspirating Accessory to the Concentrating Device. Slowly aspirate the BMA plasma into the 60mL syringe. Aspirate until the pipe completely fills with BMA RBC. Open the stopcock to the 3mL syringe and aspirate 2mL of BMA RBCs. Transfer the BMA plasma and RBC mixture into the Concentrating Accessory. Transfer the BMA RBC first.
7. **SECOND SPIN:** Counterbalance the Concentrating Accessory and place back into the centrifuge. Close centrifuge lid. Set RPM to 3800 and the time to 7 minutes. Press the start button. Once the centrifuge stops, remove the Concentrating Accessory.
8. **SECOND EXTRACTION:** Using a 60mL syringe, aspirate the BMA plasma leaving 7mL of PureBMC® in the Concentrating Accessory.
9. **RESUSPEND THE BMC:** Attach the 12mL syringe and gently swirl the Concentrating Accessory to re-suspend the bone marrow cell concentrate into the plasma.
10. **EXTRACT BMC:** Tilt the Concentrating Accessory and aspirate the PureBMC® into the 12mL syringe. Remove sterile syringe and apply sterile cap.

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

Single  
Use Only

Store between  
15 °C and 35 °C

0653

Authorized  
Representative  
Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

Australian Sponsor  
Emergo Australia  
Level 20  
Tower II, Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

EmCyte Corporation.  
4331 Veronica S. Shoemaker  
Blvd., Suite 4  
Fort Myers, FL 33916, USA  
Phone: 239-481-7725

**EmCyte Corporation**  
4331 Veronica S. Shoemaker Blvd., Suite 4  
Fort Myers, FL 33916

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






**BC60-PURE \*\* PLEASE DISCARD RED VENTED CAP FROM CONCENTRATING DEVICE BEFORE USE \*\***

**Note:** Always swab self sealing port with sterile alcohol prior to accessing with a sterile syringe.

<p>Step 1:</p>  <p>Use the filter needle to draw 10mL of Sodium Citrate Anticoagulant or Heparin 1000 units/mL into the 60mL Syringe</p>	<p>Step 2:</p>  <p>Draw 50mL of bone marrow aspirate from the patient, filling the syringe to 60mL</p>	<p>Step 3:</p>  <p>Filter from syringe to syringe by simultaneously pushing and pulling each syringe. Filter the BMA in the direction indicated on the filter.</p>	<p>Step 4:</p>  <p>Inject anticoagulated &amp; filtered BMA into the Concentrating Device</p>	<p>Step 5:</p>  <p>Counterbalance and process the Concentrating Device at:</p> <p><b>2.5 minutes</b> <b>3800 RPM</b></p> <p>Using the Executive Series Centrifuge</p>
<p>Step 5:</p>  <p>Aspirate the BMA plasma into the 60mL syringe. Then open the stopcock to the 3mL syringe and aspirate 2mL of BMA red blood cells</p>	<p>Step 6:</p>  <p>Inject 2mL BMA RBC FIRST</p> <p>Transfer the BMA RBC from the 3mL syringe first, then transfer BMA plasma from the 60mL syringe into the Concentrating Accessory</p>	<p>Step 7:</p>  <p>Counterbalance and process the Concentrating Accessory at:</p> <p><b>7 minutes</b> <b>3800 RPM</b></p> <p>Using the Executive Series Centrifuge</p>	<p>Step 8:</p>  <p>Bone marrow cell concentrate buffycoat separates out at the bottom of the Concentrating Accessory</p>	
<p>Step 9:</p>  <p>Aspirate BMA Plasma Leave approximately 7mL</p> <p>Aspirate BMA plasma from the Concentrating Accessory. Leave 7mL of plasma.</p>	<p>Step 10:</p>  <p>Attach the 12mL syringe and swirl to resuspend the BMA buffycoat into the plasma.</p>	<p>Step 11:</p>  <p>Tilt to immerse the Aspirating Pipe into the PureBMC®</p>	<p>Step 12:</p>  <p>Extract the PureBMC® into the 12mL syringe.</p>	